

# GlobalMed Camera Specification Overview

GlobalMed offers a wide selection of exam cameras and otoscopes. The following table is provided to help you better understand the unique capabilities of our most popular models.

	TotalExam Lite	TotalExam 3	Aurical OtoCam 300
SPECIFICATION			
<b>Overview</b>	High Resolution, Lightweight, Affordable. Perfect for initial entry into telemedicine or electronic records.	High-Definition imaging and versatile design for specialty use.	Affordable all-in-one design digital otoscope in a small, easy to hold design.
<b>Resolution</b>	1080 x 720	1080 x 720	720 x 720
<b>General Exam</b>	Camera Included	Camera Included	N/A
<b>Otoscope</b>	Otoscope Included	Otoscope & General Exam Heads Included	Otoscope Included
<b>Video</b>	MJPEG (30 fps)	MJPEG 60 fps)	NTSC & PAL 24 fps
<b>White Balance</b>	Auto / Manual	Auto / Manual	N/A
<b>Image Capture</b>	Intelligent Frame Capture (Best of 8 images)	Intelligent Frame Capture (Best of 16 Images)	2 Freeze Frame Buttons - Left/Right Ear
<b>Illumination</b>	4 LED / 2 LED Oto	8 LED / 2 LED Oto	2 LED
<b>USB connection</b>	USB 2.0	USB 3.0	USB 2.0
<b>Camera Status Indicator</b>	Yes	No	No
<b>Construction</b>	Drop Proof Poly Carbonite Blend	Metal / Plastic	Plastic
<b>Weight</b>	3.11 oz (88.4 g)	10 oz	8 oz / 225 g
<b>Size (Oto)</b>	190 mm	274 mm	170 mm
<b>Size (General Exam)</b>	190 mm	252 mm	N/A
<b>Regulatory</b>	EMI, EMC, CE FDA cleared and FCC approved	EMI, EMC, CE FDA cleared and FCC approved	EMI, EMC, CE FDA cleared and FCC approved
<b>Warranty</b>	1 year	1 year	1 year
<b>Optimized Color Temperature</b>	Yes	No	No
<b>Software Compatibility</b>	Windows 7, 10 & Higher MAC OS	Windows 7, 10 & Higher MAC OS	Works with eNcounter 3 / CapSure
<b>Dedicated Use Case Imager or Single imager with Oto Attachments</b>	Exam Camera + Otoscope	Exam Camera with Otoscope and General Exam head attachments	Otoscope
<b>Other</b>	Kit includes general exam camera and otoscope	Pivoting Head	Includes distal tip preheater and Cradle with magnetic switch
<b>KIT CONTENTS</b>			
<b>Polarizing Hood</b>	Included	Included	N/A
<b>Non-Polarizing Hood</b>	Included	Included	N/A
<b>One-Touch Collars</b>	Included	Included	N/A
<b>Tongue Depressor Attachment</b>	Designed into Camera	Included	N/A
<b>Tongue Depressors</b>	6	5	N/A
<b>Speculas – Adult</b>	6	5	N/A
<b>Speculas - Child</b>	6	N/A	N/A