



## LiteExam®

The LiteExam® station is designed to be the easiest and most affordable telemedicine solution that sits on your desk. LiteExam® has everything you need, nothing you don't need, to add telehealth capability to any school, office or organization. The solution sits on an existing desk or counter providing a fully-integrated telemedicine solution. It is ideal for users with limited budget that want a simple, affordable solution to provide telehealth care

### Simple & Affordable

- Designed with limited budgets in mind
- No installation required, just place it on a desk
- Easy to use software and devices

### Secure and Compliant

- Lockable physical storage for devices and accessories
- HIPAA compliant software to prevent unauthorized use and keep data secure

### “Lite” design, but Fully Capable

- Simple collection of vitals and patient info
- Easily connect to a remote practitioner
- Start basic and add devices as desired

### FEATURES

- All-in-One Touch Screen Computer
- eNCounter® Telehealth Software
- Lockable Storage for Devices and Accessories
- Accommodates TotalExam® Lite Camera and Otoscope Kit
- 4 USB Ports for Flexibility and Expansion
- Adjustable Height High Resolution PTZ\*
- Integrated Cable Management
- Stowable Keyboard Design


+1 (480) 922-0044  
globalmed.com

15023 N 73<sup>rd</sup> Street  
Scottsdale, AZ 85260  
United States





LiteExam®



### Lockable Storage Drawer



MODEL NUMBER GMR1410000

Dimensions	29" H, 26" W, 14.5" D with camera
Computer / Display	23.8" 1920x1080p HP Elite One all-in-one computer
Operating System	Windows® 10
Surface	Power-coated aluminum
Keyboard	Integrated wired keyboard and mouse
Storage	Secure device storage, slide out device drawer with keyed lock.
	Dual side storage bins
Ports	Qty 4 external USB 3.0 ports
Audio	Integrated Bang and Olufsen speaker and dual array microphone
Video Conferencing	Supports secure audio/video conferencing with PTZ camera* and eNcounter telehealth software
Camera*	PTZ Camera with 12x optical zoom
Devices	Lockable drawer accommodates TotalExam Lite General Exam Camera, TotalExam Lite Otoscope, and RNK Stethoscope.

\*optional