



## TRANSPORTABLE EXAM BACKPACK

The GlobalMed® Transportable Exam Backpack (TEB) brings a new dimension to providing telemedicine to the patient in remote environments.

This lightweight, USB-powered design provides a simple and organized solution, configured for in-home care or remote provider consultations with our suite of clinical devices including:

- TotalExam® Lite, Camera & Otoscope or TotalExam® 3.2 Exam Camera With Otoscope Attachment
- Digital Stethoscope
- TotalVitals® Health Monitor
- TotalECG® (optional add-on)
- Ultrasound (optional add-on)

Our devices allow a wide array of clinical data collection, from vitals signs including blood pressure, heart rate and temperature, to heart and lung sounds. Even in the most rural areas or isolated locations, a clinician can gather video evidence for an ENT consultation, and even capture ultrasound images, expanding care to patients anytime, anywhere.



+1 (480) 922-0044  
globalmed.com

15023 N 73<sup>rd</sup> Street  
Scottsdale, AZ 85260  
United States

Shown with TotalExam 3.2 Exam Camera and Otoscope Attachment



# TRANSPORTABLE EXAM BACKPACK



Shown with TotalExam® Lite General Exam Camera and Otoscope

Items in the Transportable Exam Backpack may be ordered separately. The complete kit includes:

## CAMERA

- TotalExam® Lite General Exam Camera and Oto | Qty 1
- or
- TotalExam® 3.2 Exam Camera with Oto Attachment | Qty 1

## KIT

- Polarizing Dermatology Hood | Qty 1
- Tongue Depressor Adapter | Qty 1
- USB 2.0 Cable, 2.5 Meters | Qty 2
- TotalVitals® | Qty 1
- Digital Stethoscope | Qty 1
- TotalECG® (optional add-on) | Qty 1
- Blood Pressure Cuffs | Qty 3
- Ultrasound (optional add-on) | Qty 1
- Headphones | Qty 1
- One Touch Collars | Qty 5
- 6" Adult Tongue Depressor | Qty 5
- Adult Specula, 4.25mm | Qty 5
- Child Specula, 3.0mm | Qty 5

## OPTIONS

- Optional Add-ons*
- TotalECG® | Qty 1
- Ultrasound | Qty 1

## Part Numbers:

- 1TEB002 with GMD5505HD20 TotalExam® Lite General Exam Camera and Oto
- 1TEB002 with GMD5504HD58 TotalExam® 3.2 Exam Camera / Otoscope

## FEATURES

- Lightweight and compact bag design
- Integrated 4 port USB hub
- Dedicated laptop pocket up to 15.6"
- Custom foam insert to support peripheral devices.
- Additional storage pockets for peripheral accessories