



ClearCam® CC1200TW VIDEO CONFERENCE CAMERA

The superior optical zoom lens, and exposure mode delivers a clear image even in low light or the extreme contrast of brightness and darkness in a conference room. The camera covers a wide shooting angle and achieves high-speed, quiet, and precise positioning with smooth PTZ operations. This enterprise-class USB 3.0 full HD PTZ video conferencing camera has positioning for middle and large clinical settings. Perfect for patient to provider video communication.

Features

- Equipped with a Sony 1/2.8 inch 2MP CMOS image sensor. It delivers a maximum resolution of 1080p 60 fps. The camera features smooth, quiet and accurate robotics even at highspeed. Offering simultaneous USB 3.0, HDMI, 3G-SDI, and Ethernet video outputs for connection flexibility.
- 12x optical zoom lens, with 72.5 degree horizontal viewing angle and 40.8 degrees of vertical viewing angle along with quiet, fast, and precise positioning movements to capture any image during a conference.
- USB3.0 Output, Plug-and-Play, compatible with your preferred video conferencing software and platform including eNcounter®.
- FECC (far-end camera control) supported.
- Standard VISCA and UVC1 control protocol, PTZ control and firmware upgrade via RS232 and USB port.
- 256 preset settings for pan, tilt, zoom, focus, exposure, and white balance modes.
- When ending a video conference or patient encounter, the camera automatically enters Privacy Mode - the tilt down of the camera to hide the lens.
- Supports multiple mounting options; available to install on desk, monitor, tripods or upside down in-ceiling.



+1 (480) 922-0044
globalmed.com

15023 N 73rd Street
Scottsdale, AZ 85260
United States

ClearCam® CC1200TW VIDEO CONFERENCE CAMERA



Package Contents	Camera, Power Adapter, Remote Control, USB 3.0 cable, HDMI Cable
Sensor	1/2.8" inch high quality, CMOS sensor, 2.41 MP
Lens	12x optical zoom, 12x digital zoom (HDMI)
View Angle	72° (Horiz), 43°(Vert), 80°(Diag)
Video System	HDMI 1080p/720p 50/60 fps USB 1080p/30 fps
Focal Length	f=3.9mm(wide)~46.8mm(tele)
S/N Rate	>>50dB
Minimum Illumination	1.0 lux (F1.6, 50IRE 30 fps)
Focus	Auto/ Manual
Gain Control	Auto/ Manual
White Balance	Auto/ Manual
Exposure Control	Auto/ Manual
Shutter Speed	1/1 ~ 1/10,000 sec
Pan Rotation Angle	+170°~170°
Tilt Rotation Angle	+90°~30°
Pan Rotation Speed	0.2°~80°/sec
Tilt Rotation Speed	0.2°~80°/sec
Preset	256
Video	USB 3.0 (MJPEG/YUV2 format)
Control	RS-232 / RS-422 / Ethernet (VISCA / PELCO D)
Power Consumption	8W @12VDC
Camera Dimension	7.48" x 7.28" x 5.12" (190 x185 x138 mm)
Weight	3.1 lbs. (1.4 kg)
Color	White
	TAA-Compliant

