



RNK STETHOSCOPE

The only telemedicine stethoscope of its kind, RNK Stethoscope is a breakthrough in remote auscultation. Not all telemedicine products are designed with true telemedicine in mind.

The quality of a traditional stethoscope, with state-of-the-art technology:

- Plugs direct into any PC computer and most tablets.
- Uses patented piezo technology for picking up body sounds with least artifact.
- No requirements for having a second doctor on transmit end. Patient can use the stethoscope on themselves.
- Four unique frequency filters to help the practitioner make a proper diagnosis.
- No buttons or dials. The easiest to use electronic stethoscope available.
- No batteries needed.
- Sleek “feels-good-in-your-hand” design



+1 (480) 922-0044
globalmed.com

15023 N 73rd Street
Scottsdale, AZ 85260
United States

MODEL NUMBER GMR11100082