



CLINICALACCESS® Station

ClinicalAccess® Station boasts an innovative modular design configurable to accommodate a wide array of examination capabilities. It scales to fit your program's evolving telemedicine needs, so you can easily expand the scope of services offered. Streamline your workflows with integrated seamless video conferencing and leading-edge medical devices. The powerful eNcounter® software from GlobalMed delivers timely, secure, robust telehealth solutions in clinical settings.

Efficient and Intelligent

- Simple to launch and use
- Customizable to your workflow
- Integrated software and hardware

Secure and Compliant

- Integrated security features
- Locking compartment option, prevents unauthorized usage

Sleek and Mobile

- Sleek, ovular design
- Small footprint and five-wheel base with 5" casters for stability, ease of egress, and turning
- Adjustable ergonomics; worksurface movement is 33-44"


+1 (480) 922-0044
globalmed.com

15023 N 73rd Street
Scottsdale, AZ 85260
United States





CLINICALACCESS® Station



Dimensions	71" minimum – 77" maximum H x 27"L x 24"W; 180 lbs. minimum
AC Power	AC, 90V to 220V, 50Hz to 60Hz
Battery Power	55AH, 1 hr, 1.60V/cell, 25°C/77°F
Computer	Integrated Dell® computer or similar
Operating System	Windows® 7, 10
Connectivity	Ethernet, WiFi
Outlets	Ethernet and AC power, located on back of integrated dust cover
Ports	4 USB 3.0 ports, 1 HDMI, two 5 Volt Aux Connection, 1 Ethernet
Cables	15' Ethernet cable reel 10' Power cable reel
Display(s)	27" / 22" Dual Monitors, swivel up to 270° 27" Single Curved Monitor, swivels up to 90° 34" Single Monitor, swivels up to 90° DC powered, 1080p Adjustable height: 36" min - 42" max
Camera	PTZ camera, standard 12x optical zoom
Table Top	Dimensions: 24"w x 28"d Adjustable height 36" min – 42" max
Audio	Dell® Sound Bar 17.4"W x 2.4"D x 1.4"H Speaker Type: Active 1x sound bar 5 Watt - 90-20000 Hz
Portability	5" dual caster (5), locking caster (2)
Devices	Compatible with all GlobalMed approved medical devices
Video Conferencing	Compatible with most major hardware and software codecs including: eNcounter®View, EasyShare®, Cisco®, Polycom™, Vidyo®, Avaya®, and others. Contact the GlobalMed sales team to learn more