



## CLINICALACCESS® Station Lite

ClinicalAccess® Station Lite (CAS Lite) frees you from the space concerns of larger form factors. Customize your station to accommodate your examination needs using its innovative, modular design. With its slim footprint, built-in WiFi and Ethernet, CAS Lite is the perfect mobile telemedicine solution to seamlessly fit into small clinical spaces.

### Efficient and Intelligent

- Simple to launch and use
- Customizable to your workflow
- Integrated software and hardware

### Secure and Compliant

- Integrated security features
- Locking compartment option, prevents unauthorized usage

### Sleek and Mobile

- Sleek, small form factor design
- Small footprint and five-wheel base with 5" casters for stability, ease of egress, and turning
- Highly maneuverable in confined spaces
- Spring-assisted adjustable ergonomics; worksurface movement is 33-44"


+1 (480) 922-0044  
globalmed.com

15023 N 73<sup>rd</sup> Street  
Scottsdale, AZ 85260  
United States





# CLINICALACCESS® Station Lite



<b>Dimensions</b>	20"W x 24"L x 53"H (lowest position, excluding camera)
<b>AC Power</b>	90V to 220V, 50Hz to 60Hz Voltage: Nominal 120V Input Voltage Range: 90VAC- 148VAC Frequency: 60Hz ± 3 Hz Plug Type: NEMA 5-15P
<b>Battery Power</b>	Dual battery Range of 6-8 Hours (varies by usage) Sealed Lead-Acid (SLA) VA 1000VA, Watts 600W, Size 12V/9AH Typical Recharge Time: 8 hours
<b>Outlets</b>	Ethernet and AC power, located on back of integrated dust cover
<b>Computer</b>	Integrated Dell® or like model
<b>Ports</b>	USB hub (4 ports) 5-port switch option within base
<b>Video Conferencing</b>	Compatible with most major hardware and software codecs including: eNcounter® View, Cisco®, Vidyo®, Avaya®, and others. Contact the GlobalMed sales team to learn more.
<b>Display</b>	27" LG Monitor 22" LG Monitor Dual 22"/ 16" Monitor with swing arm
<b>Camera</b>	PTZ camera, standard 12x optical zoom
<b>Sound</b>	Dell® Sound Bar 17.4"W x 2.4"D x 1.4"H Speaker Type: Active 1 x sound bar 5 Watt - 90 - 20000 Hz
<b>Table Top</b>	19 1/2"W x 20" L; Table top height range is 33" to 44"
<b>Security</b>	Bin storage in tabletop, with combination lock
<b>Cables</b>	10' Ethernet cable reel, 8' power cable reel

MODEL NUMBER GMR14112195