



## CLINICALACCESS® Station Lite

ClinicalAccess® Station Lite (CAS Lite™) frees you from the space concerns of larger form factors. Customize your station to accommodate your examination needs using its innovative, modular design. With its slim footprint, built-in WiFi and Ethernet, CAS Lite is the perfect mobile telemedicine solution to seamlessly fit into small clinical spaces.

### Efficient and Intelligent

- Simple to launch and use
- Customizable to your workflow
- Integrated software and hardware

### Secure and Compliant

- Integrated security features
- Locking compartment option, prevents unauthorized usage

### Sleek and Mobile

- Sleek, small form factor design
- Small footprint and five-wheel base with 5" casters for stability, ease of egress, and turning
- Highly maneuverable in confined spaces
- Spring-assisted adjustable ergonomics; worksurface movement is 33-44"


+1 (480) 922-0044  
globalmed.com

15023 N 73<sup>rd</sup> Street  
Scottsdale, AZ 85260  
United States





# CLINICALACCESS® Station Lite



Dimensions	20"W x 24"L x 53"H (lowest position, excluding camera)
AC Power	90V to 220V, 50Hz to 60Hz Voltage: Nominal 120V Input Voltage Range: 90VAC- 148VAC Frequency: 60Hz ± 3 Hz Plug Type: NEMA 5-15P
Battery Power	Dual battery Range of 6-8 Hours (varies by usage) Sealed Lead-Acid (SLA) VA 1000VA, Watts 600W, Size 12V/9AH Typical Recharge Time: 8 hours
Outlets	Ethernet and AC power, located on back of integrated dust cover
Computer	Integrated Dell® or like model
Ports	USB hub (4 ports) 5-port switch option within base
Video Conferencing	Compatible with most major hardware and software codecs including: eNcounter®View, Cisco®, Vidyo®, Avaya®, and others. Contact the GlobalMed sales team to learn more.
Display	27" LG Monitor 22" LG Monitor Dual 22"/ 16" Monitor with swing arm
Camera	PTZ camera, standard 12x optical zoom
Sound	Dell® Sound Bar 17.4"W x 2.4"D x 1.4"H Speaker Type: Active 1 x sound bar 5 Watt - 90 - 20000 Hz
Table Top	19 1/2"W x 20" L; Table top height range is 33" to 44"
Security	Bin storage in tabletop, with combination lock
Cables	10' Ethernet cable reel, 8' power cable reel

